



Self-contained for steam sterilization

Product with minimum request



Description

STERILtest GST E6 self-contained

Devices with nutrient medium, pH indicator and discs with spores of *G. stearothermophilus* ($1-5 \times 10^6$ CFU/disc).

Packaging
100 units

Ref.
91100

STERILtest GST E5 self-contained

Devices with nutrient medium, pH indicator and discs with spores of *G. stearothermophilus* ($1-5 \times 10^5$ CFU/disc).

100 units

91101 #

Ampoules for steam sterilization



Description

STERIL CONTROL GST E6 ampoules

Ampoules with nutrient medium, pH indicator and spores of *G. stearothermophilus* ($1-5 \times 10^6$ CFU/ampoule).



Packaging
10 Ampoules
50 Ampoules

Ref.
91010
91050

STERIL CONTROL GST E5 ampoules

Ampoules with nutrient medium, pH indicator and spores of *G. stearothermophilus* ($1-5 \times 10^5$ CFU/ampoule).

10 Ampoules
50 Ampoules

91051
91052

Strips for steam sterilization

Product with minimum request



Description

STRIP CONTROL GST E6

Strips with spores of *G. stearothermophilus* ($1-5 \times 10^6$ CFU/strip) and tubes of nutrient broth with pH indicator.

Packaging
20 strips
20 tubes
100 strips

Ref.
91055
91155

STRIP CONTROL GST E5

Strips with spores of *G. stearothermophilus* ($1-5 \times 10^5$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips
20 tubes

91056 #

STRIP CONTROL GST E4

Strips with spores of *G. stearothermophilus* ($1-5 \times 10^4$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips
20 tubes

91057 #

STRIP CONTROL GST E3

Strips with spores of *G. stearothermophilus* ($1-5 \times 10^3$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips
20 tubes

91058 #



Spores suspensions for steam sterilization

Description	Packaging	Ref.
SUSPENSION SPORE GST E8 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82002
SUSPENSION SPORE GST E7 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82003
SUSPENSION SPORE GST E6 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82004
SUSPENSION SPORE GST E5 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82005

Strips for dry heat or ethylene oxide sterilization

Description	Packaging	Ref.
STRIP CONTROL BAT E8 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁸ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91059 #
STRIP CONTROL BAT E7 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁷ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91062 #
STRIP CONTROL BAT E6 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91063
STRIP CONTROL BAT E5 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	100 strips 20 strips 20 tubes	91067 91064 #
STRIP CONTROL BAT E4 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91065 #
STRIP CONTROL BAT E3 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91066 #

Spore suspension for dry heat or ethylene oxide sterilization



Description	Packaging	Ref.
SUSPENSION SPORE BAT E10 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ¹⁰ CFU/mL).	10 mL vial	82006
SUSPENSION SPORE BAT E9 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁹ CFU/mL).	10 mL vial	82007
SUSPENSION SPORE BAT E8 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82008
SUSPENSION SPORE BAT E7 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82009
SUSPENSION SPORE BAT E6 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82010
SUSPENSION SPORE BAT E5 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82011



Self-contained for ethylene oxide sterilization

Product with minimum request



Description

STERILtest BAT E6 self-contained

Devices with nutrient medium, pH indicator and discs with spores of *B. atrophaeus* ($1-5 \times 10^6$ CFU/disc).

Packaging

100 units

Ref.

91150

STERILtest BAT E5 self-contained

Devices with nutrient medium, pH indicator and discs with spores of *B. atrophaeus* ($1-5 \times 10^5$ CFU/disc).

100 units

91151 #

Crusher for self-contained

Description

CRUSHER

Steel pincer for the easy break of self-contained vials.

Packaging

1 pincer

Ref.

96790

Strips for irradiation sterilization

Product with minimum request

Description

STRIP CONTROL PUM E7

Strips with spores of *B. pumilus* ($1-5 \times 10^7$ CFU/strip) and tubes of nutrient broth with pH indicator.

Packaging

20 strips

Ref.

91071 #

20 tubes

STRIP CONTROL PUM E6

Strips with spores of *B. pumilus* ($1-5 \times 10^6$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips

91072

20 tubes

STRIP CONTROL PUM E5

Strips with spores of *B. pumilus* ($1-5 \times 10^5$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips

91073 #

20 tubes

STRIP CONTROL PUM E4

Strips with spores of *B. pumilus* ($1-5 \times 10^4$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips

91074 #

20 tubes

STRIP CONTROL PUM E3

Strips with spores of *B. pumilus* ($1-5 \times 10^3$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips

91075 #

20 tubes

Steel coupons for hydrogen peroxide sterilization

Product with minimum request



Description

OXI CONTROL E6 SS

Steel coupons with spores of *G. stearothermophilus* ($1-5 \times 10^6$ CFU/coupon) and tubes of nutrient broth with pH indicator.

Packaging

20 coupons

20 tubes

Ref.

91082

OXI CONTROL E5 SS

Steel coupons with spores of *G. stearothermophilus* ($1-5 \times 10^5$ CFU/coupon) and tubes of nutrient broth with pH indicator.

20 coupons

20 tubes

91083 #



Poliflex coupons for hydrogen peroxide sterilization

Product with minimum request



Description

OXI CONTROL E6 PLX

Poliflex coupons with spores of *G.stearothermophilus* ($1-5 \times 10^6$ CFU/coupon) and tubes of nutrient broth with pH indicator.

Packaging

20 coupons
20 tubes

Ref.

91085

OXI CONTROL E5 PLX

Poliflex coupons with spores of *G.stearothermophilus* ($1-5 \times 10^5$ CFU/coupon) and tubes of nutrient broth with pH indicator.

20 coupons

20 tubes

91086 #

Strips for hydrogen peroxide sterilization

Description

OXI CONTROL E6 STRIP

Strips with spores of *G.stearothermophilus* ($1-5 \times 10^6$ CFU/strip) and tubes of nutrient broth with pH indicator.

Product with minimum request

Packaging

20 strips
20 tubes

Ref.

91088

OXI CONTROL E5 STRIP

Strips with spores of *G.stearothermophilus* ($1-5 \times 10^5$ CFU/strip) and tubes of nutrient broth with pH indicator.

20 strips

20 tubes

91089 #

Spore suspensions for hydrogen peroxide sterilization

Description

SUSPENSION SPORE GST VHP E8

Sterile suspension of *G. stearothermophilus* spores in distilled water ($1-5 \times 10^8$ CFU/mL).

Packaging

10 mL vial

Ref.

82030

SUSPENSION SPORE GST VHP E7

Sterile suspension of *G. stearothermophilus* spores in distilled water ($1-5 \times 10^7$ CFU/mL).

10 mL vial

82031

SUSPENSION SPORE GST VHP E6

Sterile suspension of *G. stearothermophilus* spores in distilled water ($1-5 \times 10^6$ CFU/mL).

10 mL vial

82032

SUSPENSION SPORE GST VHP E5

Sterile suspension of *G. stearothermophilus* spores in distilled water ($1-5 \times 10^5$ CFU/mL).

10 mL vial

82033

PVC coupon for hydrogen peroxide sterilization

Description

OXI CONTROL E6 PVC

PVC coupons with spores of *G.stearothermophilus* ($1-5 \times 10^6$ CFU/coupon) and tubes of nutrient broth with pH indicator.

Packaging

20 coupons
20 tubes

Ref.

91092

Aluminum coupons for hydrogen peroxide sterilization

Description

OXI CONTROL E6 ALUMINUM

Aluminum coupons with spores of *G.stearothermophilus* ($1-5 \times 10^6$ CFU/coupon) and tubes of nutrient broth with pH indicator.

Packaging

20 coupons
20 tubes

Ref.

91098

Glass coupons for hydrogen peroxide sterilization

Description

OXI CONTROL E6 GLASS

Glass coupons with spores of *G.stearothermophilus* ($1-5 \times 10^6$ CFU/coupon) and tubes of nutrient broth with pH indicator.

Packaging

20 coupons
20 tubes

Ref.

91096

